

Rhythm Multi 3-1 Alternating System

Rhythm Multi 3-1 Alternating System



Main Feature:

- * 8" air cell mattress replacement
- * 3-1 alternation and static
- * Low air loss therapy
- * Pulsation therapy
- * Happy Heel – dial softness control at foot section
- * Auto Firm provides immediate firm surface
- * High vapor permeable and quilted nylon cover
- * 2" convoluted foam base
- * Cell in cell design provides addition protection for upper torso and sacrum during power outage
- * Integrated glide sheet to base cover for easy transferring
- * Visual and audio alarm for low pressure and power failure with alarm mute function
- * Keypad lock out function
- * 10 Digital comfort setting
- * Modularized design for air mattress
- * Integrated power cable retainer
- * SWL 600lb
- * Quick CPR deflation within 20 seconds

Clinical Considerations:

- * Multiple stage 2, 3 or 4 pressure ulcer

Rhythm Multi Specification

Control Unit

Model Name	<i>Rhythm Multi</i>
Model No.	FC-PHR0009

Size (inch) LxWxH	17.7" x 6.8" x 10.8"
Weight (lbs)	13.2
Cycle Time	3~95 min
Min Operating Pressure	10 +/- 3mmHg
Max Operating Pressure	61 +/- 4mmHg
Flow-rate	N/A
Rated Voltage	AC 110-120V
Rated Frequency	60 Hz
Fuse Rating	5A 250V
Max Current	5A
Classification	Class I, Type BF Not AP or AGP type
Operation Temperature	15°C (59°F)~ 35°C(95°F)
Operation Humidity	30% ~ 75%
Mode of Operation	Continuous
Standard	IEC 60601-1, CAN/CSA C22.2 No. 601.1, IEC 60601-1-2

Mattress

Model No	FM-PHR0007
Size (inch)	36" (W) x 80" (L) x 10" (H)
Weight (lbs)	28.6
Cell Number	18 cells
Cells Material	Nylon coated with PU
Cover Type	Zipper cover
Cover Material	Nylon woven fabric w/ PU coating finish
Base Material	Woven Polyester fabric w/ PVC backing