

Rhythm Turn System

Rhythm Turn System



Main Feature:

- * 8" air cell mattress replacement
- * Low air loss therapy
- * 20 and 40°C lateral rotation with side perimeter and static therapy
- * Auto Firm provides immediate firm surface
- * High vapor permeable and quilted nylon cover
- * 2" convoluted foam base
- * Integrated glide sheet to base cover for easy transferring
- * Visual and audio alarm for low pressure and power failure with alarm mute function
- * Keypad lock out function
- * 10 Digital comfort setting
- * Modularized design for air mattress
- * Integrated power cable retainer
- * SWL 600lb
- * Quick CPR deflation within 20 seconds

Clinical Considerations:

- * Acute intervention for pulmonary patients

Rhythm Turn Specification

Control Unit

Model Name	<i>Rhythm Turn</i>
Model No.	FC-PHR0010
Size (inch) LxWxH	17.7" x 6.8" x 10.8"
Weight (lbs)	13.8

Cycle Time	3~95 min
Min Operating Pressure	10 +/- 3mmHg
Max Operating Pressure	61 +/- 4mmHg
Flow-rate	N/A
Rated Voltage	AC 110-120V
Rated Frequency	60 Hz
Fuse Rating	5A 250V
Max Current	5A
Classification	Class I, Type BF Not AP or AGP type
Operation Temperature	15°C (59°F)~ 35°C(95°F)
Operation Humidity	30% ~ 75%
Mode of Operation	Continuous
Standard	IEC 60601-1, CAN/CSA C22.2 No. 601.1, IEC 60601-1-2

Mattress

Model No	FM-PHR0013
Size (inch)	36" (W) x 80" (L) x 10" (H)
Weight (lbs)	37.4
Cell Number	18 cells
Cells Material	Nylon coated with PU
Cover Type	Zipper cover
Cover Material	Nylon woven fabric w/ PU coating finish
Base Material	Woven Polyester fabric w/ PVC backing